STARSWAB® ANAEROBIC TRANSPORT SYSTEM

INTENDED USE

Starswab® Anaerobic Transport System is a self-contained, sterile, ready-to-use system, which includes a tube containing a semi-solid, pre-reduced transport medium. It is intended for maintaining the viability of anaerobic as well as aerobic and facultative anaerobic microorganisms during transportation for up to 48 hours. Sterile packaging allows for use in sterile fields.

DESCRIPTION

Once a clinical specimen has left a body site any anaerobic microorganism present must be protected from the toxic effect of atmospheric oxygen until the specimen is properly set up anaerobically.¹ Therefore, it is wise to routinely use optimal transport procedures. The use of a “gassed” out collection vial is strongly recommended since it will provide an anaerobic environment during transportation of the specimens to the laboratory.² It is also reported that the viability of the aerobic and facultative microorganisms is enhanced in an oxygen free environment.³ Viability was maintained for at least 48 hours at 22-25°C (72-77° F).⁴

MEDIUM FORMULATION

Modified Cary Blair medium supplemented with L-cysteine and Resazurin. The medium is balanced, reduced and non-nutritive. It is devised to improve the survival of clinical pathogens during transportation.⁵

PRINCIPLE OF THE PROCEDURE

Starswab® Anaerobic Transport System is an enclosed system recommended for anaerobic transportation of clinical specimens. The medium contains Sodium Thioglycollate as a reducing agent. The Resazurin indicator in the medium will turn pink to lavender if the anaerobic condition in the tube is compromised. However, a pink to lavender band (up to 2 mm) on the top portion of the medium in an unopened tube is acceptable. The recessed butyl rubber stopper and screw cap facilitate the transfer of liquid specimens to the vial through the rubber stopper using syringes or needles.

Starswab® Anaerobic Transport System

Starswab® Anaerobic Transport Tube with Double Sterile Swabs	10/Box	10 Boxes/Case	Cat. No. S120-D (7 ml) Cat. No. S120-5 (5 ml)
--	--------	---------------	--

SPECIMEN COLLECTION AND TRANSPORTATION

Peel open sterile package. Remove swab and collect specimen. Insert swab into tube. Break the scored shaft of the swab. Replace and tighten the cap quickly. Transport to the laboratory. Liquid specimens should be injected into the tube through the rubber stopper. Information on specimen collection and transport may be found in standard reference materials.⁶

MATERIAL PROVIDED

- Glass tube containing pre-reduced, buffered, semi-solid Cary Blair medium supplemented with L-cysteine and Resazurin.
- Two sterile rayon tipped swabs with scored plastic shafts.

MATERIAL NOT PROVIDED

- Syringes and needles for liquid specimen collection.
- Forceps.
- Material for identification and for other microbiological procedures.

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WARNINGS AND PRECAUTIONS

This product is for single use only. Reuse may cause risk of infection, cross contamination and/or inaccurate results.

Use aseptic procedures throughout specimen collection, culturing and identification.

Handle biological samples and needles carefully according to safe laboratory practices. In the event of any exposure to biological samples obtain appropriate medical attention.

All specimens, tubes and media should be incinerated or autoclaved prior to disposal in accordance with state and local regulations.

Expiry date on media tube represents the expiry date of the product.

LIMITATION OF USE

Microorganisms survive only for a limited time in transport media and should be transported to the laboratory as soon as possible. Avoid extremes in temperature during transportation. Specimens with less than 100,000 CFU/ml may not survive for more than 24 hours.

QUALITY CONTROL

Each lot of Starswab® Anaerobic Transport System is subjected to the following Quality Control evaluations:

1. Sterility - Representative samples of each lot are tested for sterility using Fluid Thioglycollate Broth and Tryptic Soy Broth.
2. Performance - Samples of each lot are challenged by standard microorganism inocula and tested for recovery on appropriate plated media according to current CLSI, Quality Control of Microbiological Transport Systems, M40 - A2.
3. Physical appearance - Medium level, medium integrity, pink to lavender line formation and label information are some of the physical attributes of the products that are inspected.
4. pH - Each lot of product is checked against the standard pH of 7.25 +/- 0.25.

REFERENCES

1. Haldeman, L.V. and Moore, W.E.C., Editor: Anaerobic Laboratory Manual, ed. 3, Blacksburg, Va, 1975, Virginia Polytechnic Inst. and State Univ.
2. Attebery, H.R., and Finegold, S.M.: Combined Screw-Cap and Rubber-Stopper Closure for Hungate Tubes (Pre-Reduced, Anaerobically Sterilized Roll Tubes) Appl. Microbiology, 18: 558-561.
3. Finegold, S.M., Sutter, V.L., Attebery, H.R. and Rosenblatt, J.E., 1974. Isolation of Anaerobic Bacteria, pp. 365-375. In: E.H. Lennette, E.H. Spaulding and J.R. Tovan (ed) Manual of Clinical Microbiology. ed. 2. Am. Soci. Microbiology, Washington, D.C.
4. Sutter, V.L., Attebery, H.R., Rosenblatt, J.E., Bricknell, K.S., and S.M. Finegold. 1972. Anaerobic Bacteriology Manual. Extension Division. University of California. Los Angeles.
5. Jorgensen, Pfaller, Carroll, Funke, Landry, Richter, Warnock, Manual of Clinical Microbiology, 11th Edition, 2015.

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