MULTITRANS™ SYSTEM

INTENDED USE

The Multitrans™ System is used for the collection and transport of viruses, chlamydiae, mycoplasmas and ureaplasmas obtained from clinical specimens. It is intended as an aid in the diagnosis of diseases associated with infectious organisms.

SUMMARY AND EXPLANATION

The Multitrans™ System is a self-contained, ready to use system that allows for collection and safe transport of clinical samples from collection site to the testing laboratory. The medium is osmotically balanced and buffered to maintain viral, chlamydial, mycoplasmal, and ureaplasma culture. The addition of serum albumin provides a stabilizing effect. The high concentration of sucrose aids in the preservation of viruses and chlamydiae if the specimens are frozen (-70°C /-94°F) for prolonged storage time^{2,3}. The uniquely designed cap engages the swab and secures it upon reopening. The glass beads aid in the disruption of cell associated viruses, releasing the viruses into the surrounding medium, and subsequently increasing the potential for viral recovery. The plastic swabs have either polyester or nylon flocked tips.

STORAGE OF KIT

Optimum storage temperature range is 4°C to 25°C (40°F to 77°F).^{5,6}

WARNINGS AND PRECAUTIONS

1. This product is for single use only. Reuse may cause risk of infection, cross contamination and/or inaccurate results
2. For in-vitro diagnostic use only (in the U.S.A)
3. Delays in transportation may reduce recovery of organism
4. Do not use when transport medium changes color or appears turbid.
5. Do not use after expiry date.
6. Use aseptic technique when collecting and handling specimen.
7. Observe universal precautions
8. There is always a potential for the presence of blood borne viruses including human immunodeficiency virus and hepatitis virus in the specimens. Special precautions should be taken when handling specimens that may have come in contact with blood and other bodily fluids. Follow state, local and institutional guidelines for the handling and disposition of this and all biohazard waste.

MATERIAL PROVIDED

This transport system contains

1. Plastic vial containing medium and glass beads.
2. Scored sterile swab applicator with polyester tip or nylon flocked tip.
(Visit our web site www.starplexscientific.com for swab type and packaging format.)

MATERIAL NOT PROVIDED

Other materials required for culturing, identification and titering purposes.

SPECIMEN COLLECTED

Proper specimen collection from patient is extremely critical for successful isolation and identification of infectious organisms. Specimens should be collected as soon as possible after onset of disease.

1. Peel open the pouched seal pack.
2. Aseptically remove sterile swab from package.
3. Collect specimen according to protocol.
4. Aseptically remove cap from vial.
5. Insert swab into medium and break swab at scored line.
6. Replace cap to vial. Close tightly.
7. Fill out vial label with patient information.

SPECIMEN TRANSPORTATION AND STORAGE

Transport the specimen to the laboratory as soon as possible. Best recovery is obtained when specimens are refrigerated at 2°C to 8°C (36°F to 46°F) or kept on wet ice following collection and while in transit. If there is a long delay before processing specimens should be frozen at -70°C (-94°F). Specific requirements for the shipment and handling of specimens should be in full compliance with state, federal or local regulations.

IN HOUSE QUALITY CONTROL

Representative samples of each lot are tested for sterility using Fluid Thioglycollate Broth and Tryptic Soy Broth.

Medium is checked against the standard pH of 7.2-7.8.

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LIMITATIONS

This product is not intended to be used for the collection and transportation of bacterial cultures. Multitrans™ System is not intended to replace tissue culture medium for isolation of chlamydia and viruses. Repeated freezing and thawing of specimen may reduce the number of infectious organisms.

PERFORMANCE CHARACTERISTICS

The following tests were carried out to demonstrate that the Multitrans™ System is comparable to a predicate device in its performance characteristics.

For viral and chlamydia comparison, 1 ml aliquot of virus and chlamydia preparations was added to 10 ml of Multitrans TM medium and the predicate device. Each preparation was thoroughly mixed and split into 2 equal volumes. One set of these volumes was placed at 4°C and the second at 25°C (RT).

At 0, 24, 48 and 96 hours an aliquot was removed from each vial and titered either by the plaque assay or fluorescence focus forming assay.

For Mycoplasma and Ureaplasma, strains of the organisms were grown to simulate patient specimen titre and then inoculated into transport media and subsequently inoculated onto plates at 0, 24, 48 and 72 hours to determine their viability.

Clinical study results will be available upon request.

REFERENCES

- George E. Kenny, Mycoplasma, p. 407-411, Manual of Clinical Microbiology, 4th ed., American Society for Microbiology, 1985.
- Mahony, J.B. and Chernesky, M.A. 1985, J. Clin. Microbiol. 22: 865 - 867.
- Murray, Baron, Pfaller, Tenover, Tenover, Tenover, Yolken. Manual of Clinical Microbiology, 6th ed., American Society for Microbiology, 1995.
- Blackwell, G. Jamieson, F. Fearon, M. Burton, L. Rohoman, K. Bosy R. Evaluation of Starplex Multitrans Transport Media S160 Stored at Room Temperature vs 4°C. Ontario Ministry of Health and Long Term Care, Laboratories Branch. 2002.
- Blackwell, G. Jamieson, F. Fearon, M. Burton, L. Rohoman, K. Bosy R. Bruenza, G. Evaluation of Three Commercial Transport Media Compared to In-house Transport Media, Demonstrated by the Recovery of Viruses, Chlamydiae and Mollicutes.
- R. Wheeler, D. Hathaway, T. Hayes, D. Jaskot, J. Krahn, R. Lannigan, A. Tirolese. Importance of the Starplex S160-M (Swab and Viral Transport) in the Diagnosis of the 2003 - 2004 Influenza A Outbreak. London Laboratory Services Group, London, On. Canada.
- R. Wheeler, D. Hathaway, T. Hayes, D. Jaskot, J. Krahn, A. Tirolese, R. Lannigan. The Evaluation of Starplex S160 Viral Transport Medium for Long Term Storage of Patient Specimens for Research and Development Activities in a Clinical Diagnostic Virology Laboratory. London Laboratory Services Group, London, On. Canada.
- J. Harrison, M. Garassi, S. Wicrzicki, W. LeBar. Use of the Starswab Multitrans Collection and Transport System for Detection of *Chlamydia trachomatis* (CT) and *Neisseria gonorrhoeae* (GC) with APTIMA® COMBO 2 TM, APTIMA CT and APTIMA GC Assays. Hospital Consolidated Laboratories-Providence Hospital, Southfield, ML USA.
- M. Chernesky, S. Castriciano, M. Smieja, D. Jang. A Universal Transport Medium and Flocked Swabs Enhanced the Detection of *Chlamydia trachomatis* and *Neisseria gonorrhoeae* in Nucleic Acid Amplification Assays. McMaster University, St. Joseph's Healthcare, Hamilton, On. Canada.
- Anthony Tirolese. The Evaluation of Starplex S160-M Viral Transport Medium for Long Term Storage of Patient Specimens for Research and Development in a Clinical Diagnostics Virology Laboratory. London Laboratory Services Group, London, Ontario, Canada.
- D. Hathaway, T. Hayes, D. Jaskot, J. Krahn, A. Tirolese, and R. Lannigan. Evaluation of the Starplex Multitrans Media For Maintaining The Viability of Clinical Samples. Virology Department, London Laboratory Services Group, London, Ontario, Canada.
- D. Jaskot, T. Hayes, D. Hathaway, J. Krahn, S. Milburn, A. Tirolese, R. Wheeler, and R. Lannigan. Comparison of the Standard Direct Fluorescent Antibody Testing and Viral Culture to ID-Tag™ RVP from Luminox® Molecular Diagnostics. Virology Department, London Laboratory Services Group, London, Ontario, Canada.
- R. Wheeler*, D. Hathaway, T. Hayes, D. Jaskot, A. Tirolese, and R. Lannigan. Comparison of the Starplex Flocked Swab or Dacron Swab for Epithelial Cell Recovery from Nasopharyngeal Samples in Healthy Volunteers. Virology Department, London Laboratory Services Group, London, Ontario, Canada.
- R. Wheeler*, D. Hathaway, T. Hayes, D. Jaskot, J. Krahn, S. Milburn, A. Tirolese, and R. Lannigan. Integration of Rapid Cell Culture and Molecular Analysis in a Clinical Virology Laboratory. Virology Department, London Laboratory Services Group, London, Ontario, Canada.
- L.Sullivan¹, C.Kern¹, R.Milk, and T.Karnauchow^{1,2}. Rapid Detection of Influenza by Binax Now® Influenza A&B and Meridian Tru Flu Assays in a Mixed Patient Population ¹Eastern Ontario Regional Virology Laboratory, Children's Hospital of Eastern Ontario and ²University of Ottawa, Ottawa, ON, Canada.
- L.Sullivan¹, C.Kern¹, R.Milk, and T.Karnauchow^{1,2}. Performance of Meridian TruTM and Binax Now® Assays for the Detection of Rsv in Pediatric and Adult Patients ¹Eastern Ontario Regional Virology Laboratory, Children's Hospital of Eastern Ontario and ²University of Ottawa, Ottawa, ON, Canada.
- D.Hathaway, T. Hayes, D. Jaskot, J. Krahn, S. Milburn, A. Tirolese, R. Wheeler, and R.Lannigan. Evaluation of Rapid Identification Assays for Viral Markers and Determining Optimal Sample Submission Virology Department, London Laboratory Services Group, London, Ontario.

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