STARSWAB II® BACTERIOLOGY CULTURE COLLECTION AND TRANSPORT SYSTEM PACKAGE INSERT – 50 SWABS

PROTECT FROM LIGHT
Liquid Amies®
Liquid Stuart

DO NOT FREEZE
Amies Gel

STERILE R

Note: Sodium Thioglycollate

Sodium Thioglycollate is a chemical component of Amies gel, Liquid Amies, and Liquid Stuart Medium and plays an important role in maintaining organism viability. Sodium Thioglycollate has a natural sulfur-like odour that may be detected when first opening the swab pouch. This is a normal characteristic of the product.

STARSWAB II® CULTURE COLLECTION AND TRANSPORT SYSTEM

INTENDED USE

The Starswab II® Culture Collection and Transport Systems are intended for the collection and transportation of biological specimens from the patient to laboratory where direct microscopic examination and further culturing and diagnoses can be performed.

SUMMARY OF PRODUCT, PRINCIPLE AND PROCEDURE

The Starswab II® Culture Collection and Transport Systems are sterile, ready to use units for the collection of clinical specimens. Each Starswab is comprised of a sterile peel pouch containing one or two rayon swabs on a skirted plug and a plugged polypropylene tube containing one of the following transport media: Amies gel without charcoal (SP130, Amies Clear series), Amies gel with charcoal, (SP131, Amies Charcoal series), Amies gel with reduced charcoal (SP140, Amies Plus series), Liquid Stuart in polyurethane foam (SP132, Stuart series) or Liquid Amies in polyurethane foam (SP130-L, Liquid Amies series). The swabs are available in regular plastic single, regular plastic double, minitip wire, or minitip flexible wire shaft format. Liquid Stuart Medium is also available in screw capped 10 ml tubes without foam and accompanied with flocked swabs (SP132-FL series).

The Starswab II® Culture Collection and Transport Systems are an essential component for accurate diagnosis of clinical specimens. All five media employed in The Starswab II® are non-nutritive and highly reductive². This maintains the viability of microorganisms within the specimen during transportation. All five media are recommended for collection and transportation of aerobic and facultative microorganisms. The moist and deep gel in Amies Clear, Amies Charcoal and Amies Plus provide an ideal environment for the survival of anaerobic organisms³. The incorporation of charcoal in Amies Charcoal helps to neutralize fatty acids present in the specimen that may have inhibitory effects on fastidious organisms such as *Neisseria gonorrhoeae*⁴. The Amies Clear, Liquid Stuart and Liquid Amies media can be used when direct microscopic examination is desired. Amies Plus medium allows easier microscopic examination while maintaining the neutralizing effect of the charcoal. Starswab II® with flexible wire swab is suitable for specimen collection from nasopharyngeal areas.

Specimens collected with Starswab II® Culture Collection and Transportation Systems are to be transported to the laboratory where they can be inoculated onto appropriate culture media for further testing and examination. Specimens should be transported to laboratory as soon as possible and no longer than 48 hours. Fastidious microorganisms such as *N. gonorrhoeae* should be transported within 24 hours.

MATERIAL PROVIDED

1. One sterile capped tube containing one of the four transport media.
2. Either one or two of sterile rayon tipped plastic swabs, or one mini-tipped wire swab or mini-tipped flexible wire swab fitted onto biohazard safety plug, or flocked swab in separate sterile package (SP132-FL series only).

MATERIAL NOT PROVIDED

Material for culturing, identification and other microbiological procedures.

WARNINGS AND PRECAUTIONS

1. For in-vitro diagnostic use only. (In the U.S.A.)
2. Sterility of unit assured only in unopened, intact, sealed pouches.
3. Check expiry date before use.
4. Autoclave all collection materials before discarding.
5. This product is for single use only. Reuse may cause risk of infection, cross contamination and/or inaccurate result

STORAGE

1. Optimum storage temperature range is 4°C to 25°C (40°F to 77°F).
2. Do not freeze Amies Clear gel and Amies Charcoal gel.
3. Protect from light. Store Liquid Amies and Liquid Stuart medium in foil bag to protect the polyurethane foam from discoloration. In the event the polyurethane foam has discolored due to light exposure, product performance is not affected.

SPECIMEN COLLECTION

1. Peel open sterile pack.
2. Remove plug / cap from tube.
3. Remove swab and take specimen.
4. Insert swab into tube and push in plug. For flocked swabs with 10 ml tubes, insert swab into tube, break swab at scored point and screw on cap
5. Record patient's name and information on label.
6. Send specimen to laboratory for immediate processing.

SPECIMEN TRANSPORTATION

Transport specimen to laboratory as soon as possible. Avoid extreme temperature during transportation.

QUALITY CONTROL

Each lot of Starswab II® is subjected to the following quality control evaluations:

1. Performance – Samples of each lot are challenged by standard microorganism inocula and tested for recovery on appropriate plated media according to current CLSI. *Quality Control of Microbiological Transport Systems, M40-A2*.
2. Background organisms – Microbiological staining and microscopic examination are performed on random samples. This is to ensure that there is no significant background debris or cell structure present.
3. Physical appearance – Medium level, medium integrity, pouch seal, tube seal and label information are some of the physical attributes of the product that are inspected.
4. pH – Each lot of product is checked against the standard pH of 7.25 +/- 0.25.

Sterility processing - This product has been terminally sterilized by gamma irradiation exposure and validated to achieve a Sterility Assurance Level of 10⁻⁶.

LIMITATIONS

Optimum recovery and results are achieved by direct specimen plating and/or smear preparation at the time of specimen collection.

Microorganisms survive only for a limited time in transport media and should be transported to the laboratory as soon as possible. Avoid extreme temperature during transportation. Transport times for fastidious organisms such as *N. gonorrhoeae* should not exceed 24 hours, as significant loss of viability can occur. Chlamydia, mycoplasmas and viruses require special transport systems.

REFERENCES

1. Amies C.R. A Modified Formula for the Preparation of Stuart's Medium. Canadian Journal of Public Health, July 1967, Vol. 58, 296-300.
2. Isenberg H.D., Schoenenkcht F. and Von Graeventz A. Cumetech 9, Collection and Processing of Bacteriological Specimens. Coordinating editor, S.J. Rubin. American Society of Microbiology, Washington, D.C. 1979.
3. Murray P.R. et. al. Manual of Clinical Microbiology, current edition. American Society of Microbiology, Washington, D.C. 2003.
4. Quality Control of Microbiological Transport Systems; Approved Standard. CLSI Document M40-A2.
5. David Goldfarb^{1,2}, Robert Slinger^{1,2}, Ron K. Tam³, Nicholas Barrowman⁴, and Francis Chan⁵ Journal of Clinical Microbiology, Assessment of Flocked Swabs for Use in Identification of Streptococcal Pharyngitis, Sept 2009, p. 3029-3030

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RÉFÉRENCES

1. Amies C.R. A Modified Formula for the Preparation of Stuart's Medium. Canadian Journal of Public Health. July 1967. Vol. 58, 296-300.
2. Isenberg H.D., Schoenkecht F. and Von Graeventz A. Cumetech 9, Collection and Processing of Bacteriological Specimens. Coordinating editor, S.J. Rubin. American Society of Microbiology, Washington, D.C. 1979.
3. Murray P.R. et coll. Manual of Clinical Microbiology, édition actuelle. American Society of Microbiology, Washington, D.C. 2003
4. Quality Control of Microbiological Transport Systems; Approved Standard. CLSI Document M40-A2.
5. David Goldfarb^{1,2}, Robert Slinger^{1,2}, Ron K. Tam³, Nicholas Barrowman⁴, and Francis Chan² Journal of Clinical Microbiology, Assessment of Flocked Swabs for Use in Identification of Streptococcal Pharyngitis, Sept 2009, p. 3029-3030

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1. Amies C.R. A Modified Formula for the Preparation of Stuart's Medium. Canadian Journal of Public Health. July 1967. Vol. 58, 296-300.
2. Isenberg H.D., Schoenkecht F. and Von Graeventz A. Cumetech 9, Collection and Processing of Bacteriological Specimens. Coordinating editor, S.J. Rubin. American Society of Microbiology, Washington, D.C. 1979.
3. Murray P.R. et al. Manual of Clinical Microbiology, edición actual. American Society of Microbiology, Washington, D.C. 2003
4. Quality Control of Microbiological Transport Systems; Approved Standard. CLSI Document M40-A2.
5. David Goldfarb^{1,2}, Robert Slinger^{1,2}, Ron K. Tam³, Nicholas Barrowman⁴, and Francis Chan² Journal of Clinical Microbiology, Assessment of Flocked Swabs for Use in Identification of Streptococcal Pharyngitis, Sept 2009, p. 3029-3030

Pretium Canada Packaging ULC
50A Steinway Blvd., Etobicoke, Ontario, Canada M9W 6Y3
Tel / Tél. (416) 674-7474 • (800) 665-0954
Fax / Téléc. : (416) 674-6067

Item No. 02162
N° d'article. 02162
Artículo No. 02162

Rev. 02162-1
12/2019