

VersaTube™

Just as the name implies –
this container is VERSATILE.

A unique, patented tube with a hinged cap,
VersaTube™ is ideal for metals digestion,
specimen collection, transport and storage.



Cap closes securely
and can be locked
down with the snap tab

Etched flat cap surface
provides additional space
for labeling identification

Cap can be positioned
for venting or
easily removed



TECHNICAL SPECIFICATIONS

VT 502 500/case
VT-140 140/case

Capacity 50 ml Height 4.4" Width 1.145" Cap Diameter 1.36"

METALS DIGESTION

Simplify digestion analysis,
save time and money:

- Vent samples using the hinged cap – no watch glass required
- Certified accuracy of graduations (Class A) and trace metal levels
- Certificate of Analysis in each case
- Heat resistant polypropylene homopolymer
- Compatible with most digestion blocks
- Cap closed when received – no contamination concerns

COLLECTION

From field to laboratory,
collect & transport specimens:

- Easy-flip lid for fast closure
- Clear polypropylene for sample viewing
- Etched cap surface for writing
- 5 ml – 50 ml graduations for quick measurement
- Locking tab

STORAGE

Samples, swabs, pills, spices,
store anything safely:

- Clear polypropylene for viewing
- Etched cap surface for writing
- Free-standing flat bottom for storage
- Heat resistant
- Chemical and corrosion-resistant
- No moisture absorption

Customize with colour or amber (light blocking).

VersaTube™

COLLECT > TRANSPORT > STORE

VersaTube™ is ideal for a variety of markets: from permanent sample storage to specimen collection, transport or easy access to cooking ingredients or medical supplies. . . VersaTube™ is VERSATILE.

FOOD

Fresh storage of baking ingredients and spices



MINING

From coal to gold - collect everything in between



VETERINARY

Easy access to medical supplies, pills and samples



ENVIRONMENTAL

From field to lab - collect, transport and store samples



EDUCATION

Specimens under scrutiny, stored safely



Visit our websites at www.starplexscientific.com and www.mediplast.ca to discover our wide range of products suitable for a variety of markets and applications.



Dip 'n Count® Screening Kit

Ideal for clinical or industrial applications; this screening kit features a double sided media paddle conveniently attached to a removable cap and a 30ml vial.

Mediplast Pharmaceutical Bottles

Reliable storage for dry or liquid content comes in a comprehensive range of bottle shapes, colours, caps, capacity and neck sizes.



Histology Containers

Carefully engineered for the storage and light transport of liquid specimens. Secure seals reduces leakage and evaporation.

Vials

A comprehensive range of premium quality, durable polypropylene vials that coordinate with capping options for transport and storage.




Media Tubes

Reliable, sterile collection and transport of samples ensures accurate results. Optimal design features include colour coded caps, round bottoms for easy racking and a variety of swab tips and shafts to choose from.

Specimen Containers

Complete with tamper evident caps, our innovative line of leakproof containers are ideal for the safe collection, transport and storage of virtually any liquid specimen.



info@starplexscientific.com
1-800-665-0954
www.starplexscientific.com



Starplex Scientific Inc.
Etobicoke, Ontario
Canada



Starplex Scientific Corp.
Cleveland, TN
U.S.A.