

# 30ml Collector Vials and Caps Deliver Reliable Performance.

## 30ml Fecal Collector

*Integrated fecal spork cap*

Features a combination spoon and fork device to allow easy collection of any form of fecal matter. The vial has easy to read graduations which assist in the collection of the proper sample size.

## 30ml Pinworm Collector

*Integrated Pinworm Collector Cap*

For the collection and transport of pinworm ova. This 30ml vial features a clear polystyrene paddle with light adhesive on one side, allowing the pinworm and ova to be viewed under the microscope.

## 30ml Universal Collector

*Multi-purpose, graduated, sterile*

This versatile, freestanding, skirted conical bottom vial can be used for multiple purposes including fecal, urine and sputum collection.

## 30ml Slide Holder

*Durable, leak-resistant, ideal for transport*

Ideal for transporting up to four standard size slides with or without fixative. Slides do not touch and remain protected from breakage within each groove. Punch out hole in the bottom of the vial allows venting for dry transport. This unlabeled product is suitable for use as a mailer.

## Caps

- Molded "o-ring" design to provide a secure seal
- Ridges in the cap provide an easy grip
- Integrated paddle or spork design protects against possible contamination

## Vials

- Straight sided, durable polypropylene resists cracking
- Clear, easy to read graduations
- Free standing skirted conical bottom



**STARPLEX**®

A PROUD BRAND OF PRETIUM

[www.starplexscientific.com](http://www.starplexscientific.com)

1-800-665-0954