


Trusted Solutions in Specimen Collection and Transport.

As the leader in specimen collection and transport solutions, Starplex is dedicated to the development of innovative, premium quality disposables, providing reliable specimen collection and transportation for many facets of the health care industry. Starplex strives to anticipate the evolving market requirements and work cooperatively with customers to implement unique, effective solutions.

StarCLICK 100™


100ml capacity, available in sterile and non-sterile formats. The solid, straight wall design of the container paired with the durable cap creates the optimal fit for preventing costly spills and biohazardous contamination. This container is ideal for pneumatic tube systems and can be a suitable replacement for metal capped containers in certain applications.



STERILE

PRODUCT FEATURES

- Audible double click confirms the cap is correctly applied
- Clear graduations on two sides in ounces and millilitres
- Wide mouth design
- Food grade polypropylene container and polyethylene caps
- Stackable
- Tamper evident patient label



NON-STERILE



Leakbuster™ 3

Tri-sensory alignment technology

Our Leakbuster™ 3 specimen container features a tamper evident cap that ensures proper cap application by providing visual alignment points, audible closure confirmation and distinct tactile closure confirmation. Each container is lot numbered and has clear graduations. The textured, flat cap provides additional writing surface in a stackable format. Leakbuster™ 3 is also available in a clear plastic O.R. wrap format and a Midstream Kit format with a castile wipe enclosed in a clear plastic peel pack.

STARPLEX®

A PROUD BRAND OF PRETIUM

www.starplexscientific.com
1-800-665-0954

Leakbuster™ with Molded “O-Ring” Cap

This premium quality line of containers consists of 60ml, 90ml, wide 90ml and 120ml sizes for leakproof transportation. Each polypropylene container is lot numbered and has clear graduations. Ideal for operating room use, both the 90ml and 120ml containers are packaged in a sterile, clear plastic peel pouch. The wide 90ml and 120ml containers also feature a **QuickTurn™** cap closure for easy application. The patient identification labels meet current Joint Commission standards.

Sterile containers have an orange cap and tamper tab patient label. Non-sterile containers are manufactured with yellow caps. The standard 90ml container is also available in an opaque polypropylene for discrete transport. A space saver dispenser box of 20, 90 ml containers is convenient for laboratories and Doctor's offices.



Midstream Collection

For Midstream Collection, the sterile 120ml container is conveniently packaged in a clear plastic peel pouch with a castile wipe. Also available is a sterile narrow 90ml container with funnel attached and packaged with a cap, separate patient label, and castile wipe in a Ziplock transport bag.



Magenta Lined Cap

Our original magenta, flat, lined cap is available with either a 60ml or 90ml sterile polypropylene container. The standard 90ml container is also available in an opaque polypropylene for discrete transport. The non-stick silicone layer cap liner provides additional assurance against content leakage.



www.starplexscientific.com
1-800-665-0954