



Product Line	Shelf Life
Multitrans System (S160 Series)	12 months from date of manufacture
Starswab II Bacteriology Culture Collection & Transport System with Amies Clear Gel (SP130 series) and Amies Charcoal Gel (SP131 series)	24 months from date of manufacture
Starswab II Bacteriology Culture Collection & Transport System with Liquid Amies (SP130L series) and Liquid Stuart (SP132 series)	18 months from date of manufacture
Sanicult: B30U-S	12 months from date of manufacture
Starswab Anaerobic Transport System: S120-D	18 months from date of manufacture
Dry Swabs: S09 and S09D	3 years from date of manufacture
Urine Culture & Urinalysis Tube: V102-CA	Based on expiry date of boric acid tablet specified by tablet manufacturer. Typically 3 years.
Specimen Containers with Applied Temperature Strip (TS in product code)	Based on expiry date of temperature strip specified by temperature strip manufacturer. Typically 3 years.
Specimen Containers – Sterile and Non-Sterile (no temperature strip)	No expiry date issued.